

DAC2^{MX} From Puretone

DAC2^{MX} provides an easy and effective solution to a variety of hearing losses. Due to its 2 channel, 12 band technology, and being 100% digitally programmable, DAC2^{MX} systems can be precisely tailored for an individual's hearing loss.

Available in a range of custom-made models, DAC2^{MX} technology ensures clear hearing and understanding while providing you with a hearing system that is small, comfortable and completely personal to you.

The DAC2^{MX} range also includes 3 behind-the-ear models.

Simple user controls include optional memory selector push button, program switch tones and a low battery warning indicator, meaning the hearing system can be tailored specifically for your needs. Meanwhile, standard features such as the adaptive feedback manager and microphone noise reduction will ensure DAC2^{MX} hearing systems give excellent performance when it is needed most.




Features Of DAC2^{MX}

- 100% Digital Signal Processing (DSP)
- 2 DSP channels
- Up to 12 frequency band tuning
- Up to 4 programs
- Adaptive Feedback Manager (AFM)
- Multi-Mode Compression (MMC)
- Volume control (except CIC)
- Low battery warning
- Program switch tones

(Please note: Not all features available on all models)

Certain DAC2^{MX} hearing systems are custom-made for an exact fit in your ear. There are three main types of custom device.



DAC2^{MX} Technology

The DAC2^{MX} range includes state-of-the-art technology with many unique features designed to give you the best possible hearing experience, ease of use and quality of life.



AFM - Adaptive Feedback Manager

Another advanced feature of DAC2^{MX}: the Feedback Management System is not only adaptive, but is also "Response-Secured" thus retaining all the important frequency responses for the individual hearing loss. The device can also add extra stable gain which will not distort music, multi-tonal or loud sounds.



MMC - Multi-Mode Compression

Multi-Mode technology works by constantly looking for a difference in input levels and making adjustments to the incoming sound, allowing the device to capture the complete message from even quiet sources, such as conversation with a softly spoken person.

Colour Options

Custom-Made



Rapid Fit



EnigmA BTE



PT-80 & PT-90 BTEs

