

**Aquaplus colormix**Date: 22.09.2011  
Seite 1/1

Aquaplus colormix is a two-component VPS based material, which is used for swim-plugs in direct or indirect method.

| Property     | Test Method                     | Units   | Limits          |
|--------------|---------------------------------|---------|-----------------|
| appearance   | haptic                          | -----   | paste           |
| color        | visual                          | -----   | various colours |
| setting time | haptic<br>at 35°C / 95 °F       | min:sec | max. 05:00      |
| hardness     | Shore A<br>related to DIN 53505 | Shore A | 39 - 44         |

Control tests are carried out at 23°C / 73°F room temperature and 50 % relative humidity.