

Batteries & Testers

Batteries & Testers

Hearing System Batteries
Battery Testers & Tools
General Purpose Batteries



Rayovac® Zinc Air



Mercury free. Lasting up to 30% longer than equivalent mercury free hearing aid batteries, Rayovac® Mercury Free combines market leading performance with environmentally friendly technology. Available in size 675, 13, 312 & 10. Supplied in boxes of 60 (priced per cell).



- BAT675AE-6-MF** Size 675
- BAT13AE-6-MF** Size 13
- BAT312AE-6-MF** Size 312
- BAT10AE-6-MF** Size 10

Rayovac® Zinc Air - 12 Pack



Mercury free. 4 year shelf life. Available in size 675, 13, 312 & 10. Supplied in boxes of 120 (priced per cell).



- BAT675AE-12-MF** Size 675
- BAT13AE-12-MF** Size 13
- BAT312AE-12-MF** Size 312
- BAT10AE-12-MF** Size 10

PureCell Plus - 6 Pack

Mercury free. 4 year shelf life. Available in size 675, 13, 312 & 10. Supplied in boxes of 60 (priced per cell). Manufactured by Rayovac®.



Suggested Buy



- BAT675-PLUS-MF** Size 675
- BAT13-PLUS-MF** Size 13
- BAT312-PLUS-MF** Size 312
- BAT10-PLUS-MF** Size 10

Rayovac® Cochlear



Rayovac® Cochlear Mercury Free combines market leading performance with environmentally friendly technology. Available in size 675. Supplied in boxes of 60 (priced per cell).



- BAT675CP-6-MF** Cochlear

Power One Zinc Air

Mercury free. Reliable and safe hearing aid batteries, even under adverse environmental conditions. Available in size 675, 13, 312 & 10. Supplied in boxes of 60 (priced per cell).



- BATP675-MF** Size 675
- BATP13-MF** Size 13
- BATP312-MF** Size 312
- BATP10-MF** Size 10

Power One Cochlear Zinc Air

Mercury free. Reliable and safe cochlear implant batteries, even under adverse environmental conditions. Available in size 675. Supplied in boxes of 60 (priced per cell).



- BATP675CP** Cochlear Zinc Air

Duracell® Activair® Zinc Air

Activair® Mercury free batteries offer excellent performance and are suitable for all types of hearing aids. All Activair® high power Zinc Air batteries have long-reach easy-fit tabs. Available in size 675, 13, 312 & 10. Supplied in boxes of 60 (priced per cell).



- BAT675HPX-MF** Size 675
- BAT13HPX-MF** Size 13
- BAT312HPX-MF** Size 312
- BAT10HPX-MF** Size 10

Digital Battery Tester

Digital keyring-sized battery tester with storage compartment for 2 batteries. Supplied in packs of 20.



- ACC437** 20 x Digital Battery Tester

Digital Battery Tester

Digital battery tester with LCD display. Ideally sized to fit on a keyring. Supplied individually.



- ACC227D** Digital Battery Tester

Rayovac® Magnet Stick

RAYOVAC

Allows for easy insertion and removal of hearing aid batteries. Supplied in packs of 4.



ACC406A 4 x Magnet Stick

Rayovac® Battery Caddy **RAYOVAC**

Keyring sized battery holder. Holds approximately 6 hearing aid batteries. Supplied individually.



ACC408A Battery Caddy

Wonder Tool

For removing batteries stuck in most ITE and ITC hearing aids.



ACC213 Single
ACC213-10 10 Pack

General Purpose Batteries

Panasonic

Available in AA, AAA, C, D & 9V sizes. We also supply various other button cell and watch batteries. Please contact us for more information.



BATLR6 Panasonic® AA - 4 Pack

BATLR03 Panasonic® AAA - 4 Pack

BATLR14 Panasonic® C - 2 Pack

BATLR20 Panasonic® D - 2 Pack

BAT6LR61 Panasonic® 9V - 1 Pack

Power One ACCU Plus Rechargeable Batteries



Rechargeable batteries are the most environmentally friendly products of all, due to their basic reusability. During its lifecycle each battery replaces up to 57 conventional primary batteries. Power One ACCU Plus rechargeable batteries do not contain either mercury or lead or cadmium. They are nickel-metal hydride batteries in a stainless steel housing. Their capacity remains constant even after frequent recharging with our chargers (no "memory effect").



New Product

Not every hearing aid is suited to the use of rechargeable batteries - some hearing aids have too high an energy consumption. As a rule of thumb, if your hearing aid can supply current for more than 10 days without the zinc air battery needing to be changed, it is worth considering rechargeable batteries as an alternative power source. Mercury free, reliable and safe hearing aid batteries, even under adverse environmental conditions. Available in size 675, 13 & 312.



BATP675-ACCU Size 675 - 4 Pack

BATP13-ACCU Size 13 - 2 Pack

BATP312-ACCU Size 312 - 2 Pack

Power One Pocket Battery Charger



New Product

Power One pocket charger for size 13 & size 312 ACCU Plus rechargeable batteries.



ACC476 13 & 312 Pocket Battery Charger

Power One 675 Pocket Battery Charger



New Product

Power One pocket charger for size 675 ACCU Plus rechargeable batteries.



ACC227D 675 Pocket Battery Charger