

Maintenance & Hygiene

Wax Removal Sprays

Conditioning Units

Ear Lotions

Earmould Cleaning Solutions

Brushes & Tools



PerfectDry® Lux

UV-C disinfection and drying solution for hearing systems. Features forced-air fan drying system with 360 degree UV-C lamp. Powered via mains or USB. Supplied with UK mains adaptor.

NOW DRIES IN
30 MINS



ACC430H PerfectDry® Lux
ACC430K EU Power Supply
ACC430M UK Power Supply

PerfectDry®

Electric drying solution for hearing systems. Can hold 2 BTEs with earmoulds, several ITEs, glasses etc. Features secure drying cycle with timer and automatic switch-off. Effectively dries hearing systems without condensation. Supplied with UK mains adaptor.



ACC430G PerfectDry®
ACC430J EU Power Supply

Dry & Store®

Dry & Store

A complete conditioning system for all types of hearing instruments and in-the-ear devices. While you sleep, it removes moisture, dries earwax, kills germs and deodorises your hearing instrument. Dimensions: 180 mm x 99 mm x 91 mm. Supplied with UK mains adaptor and 1 x desiccant Dry-Brik II.



ACC357A Dry & Store®
ACC357D UK Mains Adaptor
ACC357H Euro Mains Adaptor
ACC357G 3 x Replacement Dry Brik II
ACC357G-3 9 x Replacement Dry Brik II
ACC357G-10 30 x Replacement Dry Brik II
ACC357C Replacement Lamp

Dry & Store® Zephyr

Dry & Store

Miniature automatic drying and conditioning system. Removes moisture from hearing instruments overnight. Dimensions: 150 mm x 92 mm x 72 mm. Supplied with UK mains adaptor and 1 x desiccant Dry-Brik II.



ACC357J Dry & Store® Zephyr
ACC357G 3 x Replacement Dry Brik II
ACC357G-3 9 x Replacement Dry Brik II
ACC357G-10 30 x Replacement Dry Brik II

PerfectClean

The first system to clean, disinfect, and dry hearing aids with a press of a button

- Vortex cleaning action effectively removes earwax from hearing aid devices
- Helps to prevent irritation
- Maintains optimal performance of hearing aid devices
- Easy to use
- Fits all types (standard and custom-made)
- 3 actions over 90 minutes
 - Cleaning: earwax removal solution
 - Disinfecting: active anti-bacterial agent + UV-C system
 - Drying: forced-air fan drying system
- Easy to refill
- Visual indication of the solution level



PerfectClean Refill



ACC430P PerfectClean
ACC430Q PerfectClean Refill

Dry Spot Dehumidifier

Stylish container with microwaveable desiccant pillow. Designed to dry earwax and kill germs in hearing systems. Improves hearing system reliability and reduces the need for repairs.



- ACC333A** 1 x Dry Spot Dehumidifier
- ACC333A-10** 10 x Dry Spot Dehumidifier
- ACC333B** 1 x Spare Desiccant Pillow
- ACC333B-20** 20 x Spare Desiccant Pillow

Moisture Control Kit

For drying hearing systems overnight. Contains: static shield bag, desiccant and 3-spot humidity card. Available as a kit or individually.



- ACC171A** 1 x Complete Kit
- ACC171A-50** 50 x Complete Kit
- ACC171B** 1 x Static Shield Bag
- ACC171C** 1 x Desiccant
- ACC171D** 1 x 3-Spot Humidity Card

Air Puffer

For removing moisture from tubing. Available in blue.



- ACC121** 1 x Air Puffer
- ACC121-10** 10 x Air Puffer

Acu-Life® Air Puffer

Helps keep earpieces and tubing clear from moisture, dirt and debris for clear, unobstructed sound. Great for earpieces, tubes and open fit hearing aids.



- ACC435C** Air Puffer

Suction Tool

Ideal for the removal of moisture and wax from hearing instruments. Features recessed release button to prevent accidental operation. The housing separates to provide easy access for cleaning.



- ACC336** Suction Tool

PowerVac

Powerful suction tool designed to remove light debris and moisture from hearing instruments and earmolds. Supplied with 2 x AA batteries.



- ACC347A** 1 x PowerVac
- ACC347A-3** 3 x PowerVac
- ACC347B** 1 x Suction Tips
- ACC347B-5** 5 x Suction Tips

Propulse® Ear Irrigator Pack

PROPULSE®
IRRIGATOR™

Featuring a bold, modern design, the new Propulse® ear Irrigator sets new standards in ear irrigation whilst retaining the familiar and trusted features that you and your patients depend upon. The new Propulse® ear Irrigator has been designed with both the patient and the clinician in mind. The new Propulse® ear Irrigator has been ergonomically designed – for instance, the new handle design makes it easier to remove the QrX™ tips. Tried and tested key Propulse® features such as easy hands-free operation via the footswitch are maintained. A range of products & accessories that complement the Propulse® ear Irrigator are also available.

- Trusted Propulse® quality.
- Hands-free operation by use of footswitch.
- Variable pressure control for comfort and safety.
- Ergonomically designed handle.
- Easy tip removal & replacement.
- Utilises Propulse® QrX™ Technology.



ACC458A Propulse® Ear Irrigator Pack

Propulse® QRX Tips PROPULSE®

Propulse® Single Use QrX™ Tips are designed exclusively for use with the Propulse® G5 Ear Irrigator. The unique locking mechanism provides an easy fit on to the quick release handle whilst securing the tip throughout the irrigation procedure. Supplied in pack of 100.



ACC458B 100 x QRX™ Tips - White

Propulse® ProTect™ PROPULSE®

This waterproof shoulder cape helps protect the patient's clothes during the ear irrigation procedure. Supplied in pack of 36.



ACC458C ProTect™

Propulse® Mushroom Valve & Washer

The Mushroom Valve is specifically designed to prevent water flowing out of the reservoir whilst filling.



ACC458D Mushroom Valve & Washer

Propulse® Carry Case

- Soft, padded carry case designed to carry the Propulse®
- Zipped flap allows quick and easy access to the Propulse®
- Prevents damage to the Propulse® during transit
- Suitable for domiciliary visits
- Convenient shoulder strap



ACC458G Carry Case

Propulse® Cleaning Tablets

Propulse® Cleaning Tablets are specially formulated for use with Propulse® Ear Irrigators.


They ensure that the irrigator is cleaned effectively without damaging the pump and working components of the unit as some other cleaning agents can. Supplied in pack containing 200 tablets.



ACC459A Cleaning Tablets

Propulse® ProScoop

The double ended design combines a scoop for manual ear wax removal at one end, and a serrated edge at the other to allow easy attachment of cotton wool (for mopping out the ear following ear irrigation). Supplied in pack of 50.



ACC460A ProScoop

Propulse® Head Lamp

Illumination of the ear canal is essential when performing ear irrigation.

- Flexible goose neck design provides flexibility of movement
- Comfortable and light weight
- Bright white light provides excellent illumination of the external auditory meatus
- Elasticated headband fits all users
- Simple push button on/off switch



ACC461A Propulse® Head Lamp
ACC461B Propulse® Head Lamp Bulb

Propulse® Noots Ear Tank

The Propulse® reusable noots ear tank is ideal for ear irrigation and other cerumen removal procedures.



ACC462A Noots Ear Tank

Acu-Life® Hearing Aid Dehumidifier **ACU-LIFE®**

Dri-Eze hearing aid dehumidifier removes damaging daily moisture. Prolongs hearing aid life, improves sound quality and helps prevent costly repairs. Contains Dri-Eze dehumidifier and protective foam pad. No batteries or electricity needed.



ACC435A Hearing Aid Dehumidifier

Acu-Life® Wax Removal Syringe **ACU-LIFE®**

Ear wax removal 20 ml syringe for the use with water or ear wax drops. Safely clean ears with the soft tri-stream tip.



ACC435B Wax Removal Syringe

Acu-Life® Ear Irrigator **ACU-LIFE®**

The Acu-Life® ear irrigator cleaning system cleans ears and helps to prevent ear wax buildup. Tri-stream soft tip and pressure control bottle to ensure safe and effective cleaning. Bottle contains 118 ml.



ACC435D Ear Irrigator

Acu-Life® Oto-Scoop **ACU-LIFE®**

Oto-scoop with flexing tip is the safe and simple way to remove unwanted ear wax. Simply insert flexing tip into the ear and sweep the ear wax out. Safety guard protects ear drum. For adult use only.



ACC435F Oto-Scoop

Wax Remover Sprays

For effective removal of earwax from hearing instruments and in-ear monitors. Available in 50 ml and 30 ml sizes (please note: 50 ml size features compressed gas - special delivery only).



Special Delivery

50 ml
Compressed Gas
Hazardous Shipping



Special Delivery

50ml
Compressed Gas
Hazardous Shipping



30ml
Pump Action

- ACC309A** PureClean (50 ml)
- ACC309A-3** 3 x PureClean (50 ml)
- ACC309A-24** 24 x PureClean (50 ml)

- ACC309B** IEM Spray (50 ml)
- ACC448G** PureClean (30 ml)
- ACC448G-10** 10 x PureClean (30 ml)

PureClean Cleaning Wipes

Easy to use moist, hygienic cleaning wipes. Ideal for a quick cleaning or use while traveling. Supplied individually wrapped in boxes of 30 wipes.



- ACC448A** Single Pack (30 Wipes)
- ACC448A-10** 10 Packs (300 Wipes)

PureClean Comfort Gel

Earmould sealant/comfort gel for easier insertion and removal of in-ear hearing systems, earpieces and earmoulds. Supplied in 5 ml tube. Supplied individually or in a pack of 10.



- ACC448B** Single Pack
- ACC448B-10** 10 Pack

PureClean Cleaning Tablets & Cup

Used regularly they can easily remove earwax accumulated through normal wear. Tablets supplied in pack of 20 or 200. Cleaning cups supplied individually or as packs of 10.



- ACC448E** 1 x Cleaning Cup
- ACC448E-10** 10 x Cleaning Cups
- ACC448F** 20 x Cleaning Tablets
- ACC448F-10** 200 x Cleaning Tablets

PureClean Drying Capsules & Cup

Ideal for everyday use, they absorb any moisture present in hearing systems. Drying Capsules supplied in pack of 4 or 40. Drying cups supplied individually or as packs of 10.



- ACC448C** 1 x Drying Cup
- ACC448C-10** 10 x Drying Cups
- ACC448D** 4 x Drying Capsules
- ACC448D-10** 40 x Drying Capsules

PureClean Hearing Aid Starter Set

Starter set containing 6 x PureClean Cleaning Wipes, 1 x PureClean Cleaning Spray, 1 x PureClean Drying Capsule and 1 x PureClean Drying Cup. Includes zip bag, brochure and PureClean Hearing Aids sleeve.



ACC447A Hearing Aid Starter Set

PureClean Noise Protection Starter Set

Starter set containing 6 x PureClean Cleaning Wipes, 1 x PureClean Comfort Gel and 1 x Air Puffer. Includes zip bag, brochure and PureClean Noise Protection sleeve.



ACC447B Noise Protection Starter Set

Maintenance Pack

All-in-one hearing system maintenance kit. Supplied with zippa bag, dehumidifier bag, stetoclip, tubing extractor, tubing threader, blue air puffer, PureClean hygienic wipes and 5 x pre-bent tubes.



ACC302 Maintenance Pack

Medi Wipes

Alcohol based cleaning wipes for hearing systems and earmoulds. Available in a box of 100, 1000 and individually.



- ACC404-1** Individual
- ACC404** Pack of 100
- ACC404-10** Pack of 1000

Cleaning Tool Kits

Convenient pocket sized multi-tools for hearing systems. Features magnet, battery tray opener, vent cleaner, wax loop, cleaning brush, as well as 2 battery storage compartments. Available with or without digital battery tester.



- ACC398** 1 x Without Battery Tester
- ACC398-10** 10 x Without Battery Tester
- ACC421** 1 x With Battery Tester
- ACC421-5** 5 x With Battery Tester

Wonder Wipes


Anti-bacterial cleaning wipes for hands, tools and surfaces. Resistant to MRSA, Salmonella, Listeria & E. Coli. Supplied in tub of 100 wipes.



- ACC431** Single Tub
- ACC431-3** 3 x Tubs

Multi-Bristle Vent Cleaner

Vent cleaning tool with a selection of multi-sized bristles.



ACC392 Multi-Bristle Vent Cleaner

AudiSpray®

100% natural solution for ear hygiene, Audispray® gently cleanses the ears by dispersing the earwax. Audispray® for age 6 and above, and is supplied in a 50 ml bottle (about 6 months of use).



ACC424A Adult
ACC424A-PACK 9 x Adult
ACC424C Display Box

AudiClean®

AudiClean® Spray is the safer alternative to cotton buds. Made from purified seawater, AudiClean® Spray gently cleanses the ears by dispersing the earwax. Used on a regular basis AudiClean® Spray will help prevent the formation of earwax. Contents 115 ml.



ACC375E Single
ACC375E-10 Pack of 10

Auragel

Earmould sealant/lubricant for easier insertion/removal of earmoulds. Contains no oil or grease. Supplied in a 20 g tube. Available individually or in a display pack of 10.



ACC285 Single
ACC285-10 Pack of 10

Eargene®

Soothing and refreshing ear lotion. Relieves itching ears, leaving the ear cool and refreshed. Supplied in a 15 ml bottle, available individually or in a display pack of 24.



ACC286 1 x 15 ml
ACC286-24 24 x 15 ml

Earol® Olive Oil Spray

Earol® Spray has a simple yet unique delivery system capable of dispensing a discrete unit dose of olive oil into the outer ear, softening compacted earwax. Earol® can be used preventively once a week to aid ear hygiene and maintain healthy ears. Supplied in 10 ml spray.



Suggested Buy



ACC449 Single
ACC449-10 Pack of 10

£3.97
£37.68

Earol® Swim Tea Tree Oil

Earol® Swim Tea Tree Oil delivers a measured dose of a unique blend of olive oil & tea tree oil into the outer ear finely coating the auditory canal creating a water resistant barrier. Olive oil & tea tree oil has natural anti-fungal, anti-bacterial, anti-viral, and anti-infection properties and when combined with the pump action metered spray delivery system will help prevent trapped water & swimmer's ear. Supplied individually



ACC449A Earol® Swim Tea Tree Oil

Brushes & Tools



Wax Loop (ACC114)



Power Magnet &
Wax Loop (ACC117J)



Brush & Wax Loop
(ACC117C)



Brush (ACC116)



Wax Loop & Magnet
(ACC115)




CENS® ProFlex Brush & Magnet
(ACC117G)



Brush & Magnet
(ACC117A)



Brush & Child Lock
Tool (ACC116L)



Wonder Tool Battery
Remover (ACC213)



Brush & Vent Cleaner
(ACC117D)



Brush, Magnet & Wax
Loop (ACC117B)



Anticer Brush
(AN3)

- ACC114** Wax Loop
- ACC115** Wax Loop & Magnet
- ACC116** Brush
- ACC116L** Brush & Child Lock Tool
- ACC117A** Brush & Magnet
- ACC117B** Brush, Magnet & Wax Loop
- ACC117C** Brush & Wax Loop
- ACC117D** Brush & Vent Cleaner
- ACC117G** CENS® ProFlex Brush & Magnet
- ACC117J** Power Magnet & Wax Loop
- ACC213** Wonder Tool
- AN3** Anticer Brush

- ACC114-150** 150 x Wax Loops
- ACC115-100** 100 x Wax Loop & Magnet
- ACC116-100** 100 x Brushes
- ACC116L-100** 100 x Brush & Child Lock Tools
- ACC117A-100** 100 x Brush & Magnets
- ACC117B-100** 100 x Brush, Magnet & Wax Loops
- ACC117C-100** 100 x Brush & Wax Loops
- ACC117D-100** 100 x Brush & Vent Cleaners
- ACC117G-100** 100 x CENS® ProFlex Brush & Magnets
- ACC117J-50** 50 x Power Magnet & Wax Loops
- ACC213-10** 10 x Wonder Tools
- AN3-50** 50 x Anticer Brushes

PureVac

A revolutionary hearing aid vacuum system featuring a stylish aluminium enclosure, 25 bar vacuum pressure, and an all-in-one, compact design. This sleek, modern system is perfect for any hearing care professional.

Compact, modern and simple to use device for the cleaning of hearing instruments. Individual vacuum and pressure styluses allow for hygienic concentrated cleaning of a device. A timed drying chamber utilising both vacuum and conducted heat provides all over cleansing.

- 25 bar vacuum
- 50°C typical chamber heat
- 7 minute run time
- Compact design
- Easily replaceable tube
- Easily replaceable stylus
- Pressure display



Pure Vac



New Product

ACC473 PureVac

Jodi-XL

Complete clinical vacuum and desiccator system in a small, yet powerful pump. Two 2000 hour-rated vacuum pumps power this unit; one pump draws moisture and debris out through a needle vacuum and the other powers the heated desiccator. Includes a selection of different-sized needles for vents and sound bores. Powered by a 15V DC transformer, this unit can be used anywhere in the world.



ACC411A Jodi-XL

ACC411B Filters - Pack of 3

ACC411C Filters - Pack of 6